### SARG: the high resolution spectrograph of TNG

R.G. Gratton, R.U. Claudi, G. Farisato, G. Martorana, M. Rebeschini

Osservatorio Astronomico, Padova, Italy

G. Bonanno, P. Bruno, A. Calí, S. Scuderi, R. Cosentino, M.C. Timpanaro

Osservatorio Astronomico, Catania, Italy

S. Desidera

Dipartimento di Astronomia, Padova, Italy

#### Summary

The status of the high resolution spectrograph of TNG (SARG) is presented. This instrument will be in operation during early spring 2000. SARG will offer both single object and long slit (up to 30 arcsec) observing modes covering a spectral range from  $\lambda$ =0.37 up to 0.9  $\mu$ m, with resolution ranging from R=19,000 up to R=144,000. Cross dispersion is provided by means of a selection of four grisms and interference filters may be used for the long slit mode. A dioptric camera images the cross dispersed spectra onto a mosaic of two 2048 × 4096 EEV CCDs (pixel size: 13.5  $\mu$ m) allowing complete spectral coverage at all resolving power for  $\lambda < 0.8 \ \mu$ m. In order to reach a high wavelength calibration precision an iodine-absorbing cell is provided. A Distribuited Active Temperature Control System (DATCS) maintains constant the temperature of all spectrograph components at a preset value.

#### 1 Introduction

SARG is the high resolution spectrograph for the Galileo National Telescope. Istrument specifications included a high spectral resolution (maximum about  $R \sim 150,000$ ), high efficiency (peak at more than 15%), rather large spectral coverage in a single shot (two shots allowing to cover the entire accessible range from 370 to 900 nm; UV was sacrified due to to budget limitations), high stability (long term stability of 5 m/s, and possibly 1 m/s over a single night). While meeting all these requirements required a complex and sophisticated instrument, design criteria also included simplicity of use and maintenance, and respect for the TNG standards and environment. The project (approved in February 1996) was scheduled as a four year project, and within a typical budget for a 4m telescope class instrument (0.85 M Euro). To realize such a complex instrument within the allowed budget, we had to take responsibility for the optical design and integration of the instrument.

Three years after its approval, SARG is under test in its own specially devoted optical laboratory in Padova. It will be moved to La Palma at the beginning of 2000.



Figure 1: Mechanical layout of SARG. The main subsystem are labeled

SARG is a collaboration of the astronomical observatories of Padova (optics, mechanics, thermal control, integration, science verification, project management), Catania (controls, software, detectors), Trieste and Palermo (management and science verification). The iodine cell is realized as the doctoral dissertion programme of S. Desidera (Dept. of Astronomy, Un. Padova), with the collaboration of G. Favero (Dept. of Chemistry, Un. Padova) and A. Calì (Catania Observatory).

This document presents a brief description of some of the spectrograph properties and components. Details about SARG controls, software, detectors, documentation, perforamances, and the iodine cell will be given in other contribution at this meeting.

#### 2 Technical Characteristics

The general layout of SARG is shown in Figure 1. SARG is mounted on an optical table rigidly attached at the fork of TNG, about 1.4 m below the Nasmyth B focus (see Figure 1). Light reaches the spectrograph through an optical train which includes three lenses (L1, L2, and L3), and a folding mirror (FM1) which redirects light exiting horizontally along the elevation axis downward to the spectrograph location. L1 and FM1 are mechanically located within DoLoRes (the TNG low dispersion spectrograph) which permanently occupies the Nasmyth B location (see Figure 1).

In analogy with several modern high resolution spectrographs, SARG has a white pupil collimator. This design, coupled with the use of an R4 echelle, of grism cross dispersing elements, and of a large field, quite long focal length dioptric camera exploiting a large size detector, composed of a mosaic of 2  $2k \times 4k$  CCDs, pixel 13.5  $\mu$ m, allowed a very compact and simple mechanical design. Instrument size is about  $2 \times 1 \times 1$  m, very small in view of the high spectral resolution achieved by the instrument when mounted on a 3.5 m telescope. This compact design also simplifies its thermal design: a special feature of SARG is its distributed active thermal control system (DATCS), which allows the spectrograph to have a nearly constant temperature of  $20 \pm 0.5$  C.



Figure 2: SARG location on the TNG arm  ${\rm B}$ 



Figure 3: Photo of FM1 within the DoLoRes input slide

	BLUE	GREEN	YELLOW	RED
Spectral Range (nm)	369 - 518	419 - 567	462 - 795	502 - 1020
Peak efficency (%)	16	16	20	18
Min. separation (arcsec)	11.3	13.9	8.1	6.0

The preslit optics (which includes a collimated part of the optical path) and the grism cross disperser also add flexibility to SARG, allowing a variety of optical modes: multiorder short slit observations with a large spectral coverage; 30 arcsec single order long slit observations (an optical derotator allows to have an arbitrary fixed orientation of the slit projected on the sky); multiorder image slicer observations at very high spectral resolution, using a Diego modified Bowen - Walraven image slicer; high precision radial velocity observations using an I<sub>2</sub> absorbing cell.

### **3** Optics

The preslit optics (realized by Fisba Optiks) consists of three lenses: L1 transforms the f/11 TNG beam into an f/34.4 one, allowing a suitable focal position far from the elevation axis (so that SARG would not interfere with DoLoRes). A second lens (L2) located after an intermediate (free) focal position collimates the beam, and a third f/11 lens (L3) focuses the beam onto the slit (actually, a selection of 9 different apertures located on a motorized rotating table, which allows resolution from R = 17,000 up to R = 144,000). The collimated portion of the beam is used for an Abbe-Koenig inverting prism, which allows derotation over a field of about 30 arcsec (projected on sky). The preslit optics also include slides and wheels for neutral and coloured filters (allowing suppression of undesired light); the two calibrating systems of SARG (the lamp array, including three quartz iodine and a Th lamp; and the iodine absorbing cell); and a slit viewer assembly. This last have two modes of operation: direct imaging (using a Nikon objective) and wavefront analysis (using a Shack-Hartmann lenslet array), which allows to close the active optics loop.

The spectrograph shutter is located after the slit; thereafter, a second folding mirror (FM2) redirects light toward the f/11 100 mm beam white-pupil collimator (constructed by SESO). In the first pupil image created by this system, we located the 31 gr/mm R4 grating (Richardson Lab: this choice allows a quite high RS product of RS = 46,000) used in a quasi-Littrow configuration (off-plane angle of 0.7 degrees). A selection of grisms (FISBA Optik) mounted on a rotating table are located in the second pupil image, allowing a variety of cross dispersion modes (see Table 1). One further position on this table is left free, to be used for long slit observations coupled with interference filters located in the filter wheel in the preslit portion of the beam.

SARG uses a single, quite long (485 mm focal length), large field dioptric camera constructed by SESO. The detector is a mosaic of 2 EEV  $2k \times 4k$  CCDs, pixel 13.5  $\mu$ m (very similar to that in use at the OIG). Measured efficiency of this CCD as a function of wavelength is given in Figure 3



Figure 4: Efficiency of SARG CCD as a function of wavelength



Figure 5: Expected efficiency of SARG at R = 144,000 (upper panel); and the solar spectrum at this resolution, when coupled with the yellow grism



Figure 6: Layout of the SARG VME control system architecture.

### 4 Mechanics

SARG mechanics has been constructed by CINEL. It consists of three main components:

- The FM1 unit carries L1 and FM1. It is mounted on the input slide of DoLoRes
- The main spectrograph body includes the slit box (where the preslit optics is mounted), the calibration lamp unit, the slit viewer box, and the mountings for the main optical elements (collimator mirrors, echelle, grisms, camera). All these are located on a Newport optical table. The tables, and all the other elements (save for the calibration lamp and the CCD dewar) are located inside a thich thermally insulating enclosure. There are 9 motorized wheels and tables, that allow remote control of all spectrograph functions. Five of them are realized by means of commercial (OWIS) table; the remaining four are annular rotating tables realized by CINEL
- The Telescope/Spectrograph interface is a steel framework that allow proper handling and positioning of SARG at the arm of the TNG fork (fine alignment is achieved by tip/tilting a thick quartz window located within the optical path)

## 5 Controls

SARG controls are described in more details in a separate contribution at this meeting; a sketch of its overall architecture is given in Figure 7. It is a VME-system based (see Figure 6); the VME-bus directly controls the Electromare TNG standard CCD controller by means of a fiber optics link and TNG standard Atenix boards. Instrument functions are controlled by RS232 links to stand-alone controllers (three OWIS motor-controller, a National board commanding the lamps, and a Lakeshore temperature monitor). All these elements are located on a single rack, a couple of meters far from the spectrograph.



Figure 7: Layout of SARG control system



Figure 8: Run of expected temperature of SARG as a function of temperature in the Nasmyth room when the DATCS is in operation. Solid line represents average (air) temperature; dotted lines are expected minimum and maximum temperatures within the spectrograph. The dashed area represents specifications

## 6 SARG Distributed Active Thermal Control System (DATCS)

Temperature within the SARG enclosure is kept at  $20 \pm 0.5$  C by means of an active thermal control system. The design of this system was done using a computer model, tested against suitable laboratory experiments. The SARG thermal specifications were set as follows:

- Long-term stability: < 0.5 C
- Thermal inhomogeneities: < 1 C
- Rate of temperature changes: < 1 C/hr

The spectrograph passive thermal properties are:

- Weight: 750 kg
- Effective Thermal capacity: 534 KJ/C
- Heat Flux through enclosure: 10 W/C
- Thermal Constant: 15 hr
- Total heat losses: 10–200 W (depending on environment temp.)

To achieve specifications, DATCS architecture consists of two main components:

- 30 MINCO CT 198 HEATERSTAT sensorless temperature controllers coupled with thermofoil resistances
- $\bullet~8$  DT-470 silicon diode sensors controlled by a Lakeshore 208 temperature Monitor

Within this design, particulare care was devoted to DATCS calibration. To this purposes, a special Heaterstat Calibration Box was constructed, a calibration procedure was defined and strictly enforced, and a maintenance plan was defined in order to mantain conditions stable over a long period of time.

The expected performances of DATCS from our computer model are within specifications (see Figure 8). Strict procedures will be enforced to avoid that during spectrograph heating temperature gradients will not be larger than specifications.

# 7 Laboratory tests

During last autumn we had two test runs (the first in September, the second one in December). During these runs, we tested optics, mechanics, controls, software, detectors, and the thermal behavior of SARG. Note that at present only one of the two CCD's of the focal plane mosaic is available (the second one has been oredered in October, and will be available in January). Additionally, the image slicer (realized by the University College of London) is still not available, so that the maximum resolution during these test was 155,000 (nominal at order center). Spectra of calibration lamps and of the Sun (through a fiber link) were acquired. We generally find that the instrument work according specifications.

Examples of these spectra are shown:



Figure 9: Spectrum of the Th lamp, acquired with the yellow grism; the slit aperture width was 75  $\mu m,$  yielding a nominal spectral resolution of 115,000



Figure 10: Solar spectrum acquired with the yellow grism; the slit aperture width was 75  $\mu{\rm m},$  yielding a nominal spectral resolution of 115,000



Figure 11: Plot of the solar spectrum around the Na D lines, as given by the Solar Atlas by Kurucz et al. (1994: upper panel), and as obtained using SARG (lower panel) at a nominal spectral resolution of 115,000



Figure 12: The same as figure 11, but for the spectral region including the core of the Ca II K line



Figure 13: FWHM (in pixels) of the Th lines in different positions on the detector. At the nominal resolution, the projected slit width was 2.5 pixels at order centers; due to the anamorphic magnification of the echelle, the projected slit width is expected to change by about  $\pm 20\%$  at the extremes of the orders visualized on the detector

- Figure 9: portion of the spectrum of the Th lamp, acquired with the yellow grism; the slit aperture width was 75  $\mu m$ , yielding a nominal spectral resolution of 115,000
- Figure 10: portion of the solar spectrum acquired with the yellow grism; the slit aperture width was 75  $\mu$ m, yielding a nominal spectral resolution of 115,000
- Figure 11: plot of the solar spectrum around the Na D lines, as given by the Solar Atlas by Kurucz et al. (1994: upper panel), and as obtained using SARG (lower panel) at a nominal spectral resolution of 115,000
- Figure 12: The same as figure 11, but for the spectral region including the core of the Ca II K line
- Figure 13: FWHM (in pixels) of the Th lines in different positions on the detector. At the nominal resolution, the projected slit width was 2.5 pixels at order centers; due to the anamorphic magnification of the echelle, the projected slit width is expected to change by about  $\pm 20\%$  at the extremes of the orders visualized on the detector